



Recombinant Monkeypox virus Profilin (A42R)

Product Code	CSB-EP837687MHV
Abbreviation	A42R
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8V4T7
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Monkeypox virus
Purity	≥85% (SDS-PAGE)
Sequence	MAEWHKIIED ISKNNKFEDA AIVDYKTTKN VLAAIPNRTF AKINPGEVIP LITNHNILKP LIGQKFCIVY TNSLMDENTY AMELLTGYAP VSPIVIARTH TALIFLMGKP TTSRRDVYRT CRDHATRVRA TGN
Source	E.coli
Target Names	A42R
Protein Names	Recommended name: Profilin
Expression Region	1-133
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.