



Recombinant Mouse Radial spoke head 1 homolog (Rsph1)

Product Code	CSB-MP837749MO
Abbreviation	Rsph1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VIG3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSDLGSEELE EEGENDLGEY EGERNEVGER HGHGKARLPN GDTYEGSYEF GKRHGQGTYK FKNGARYTGD YVKNKKHGQG TFIYPDGSRY EGEWADDQRH GQGVYYVNN DTYTGEWFNH QRHGQGTYLY AETGSKYVGT WVHGQQEGAA ELIHLNHRYP GKFMNKNPVG PGKYVFDIGC EQHGEYRLTD TERGEEEEEE ETLVNIVPKW KALNITELAL WPTLSEEQP PPEGQQQEEP QGLTGVGDPS EDIQAEGFEG ELEPRGADED VDTRFQESQE NSYDIDQGNL NFDEEPSDLQ D
Source	Mammalian cell
Target Names	Rsph1
Protein Names	Recommended name: Radial spoke head 1 homolog Alternative name(s): Male meiotic metaphase chromosome-associated acidic protein Meichroacidin Testis-specific gene A2 protein
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.