



Recombinant Arabidopsis thaliana Uncharacterized exonuclease domain-containing protein At3g15140 (At3g15140)

Product Code	CSB-YP837802DOA
Abbreviation	At3g15140
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8W566
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASAFSAFRV SLSRISPF RD TRFSYPATLA LAHTKRIMCN SSHSVSPSPS PSDFSSSSSS SSSSPSTFSL METSENARWR PMCLYYTHGK CTKMDDPAHL EIFNHDCSKE LRVAADLER KKSQEFNFFL VIDLEGKVEI LEFPILIVDA KTMEVVDLFH RFVRPTKMSE QAINKYIEGK YGELGVDRVW HDTAIPFKQV VEEFEVWLAE HDLWDKDTDW GLNDAAFVTC GNWDIKTKIP EQCVVSNINL PPYFMEWINL KDVYLNFYGR EARGMVSMR QCGIKLMGSH HLGIDDTKNI TRVVQRMLSE GAVLKL TARR SKSNMRNVEF LFKNRIK
Source	Yeast
Target Names	At3g15140
Protein Names	Recommended name: Uncharacterized exonuclease domain-containing protein At3g15140 EC= 3.1.-.-
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.