



Recombinant *Prunus laurocerasus* Maturase K (matK)

| | |
|--------------------------|--|
| Product Code | CSB-BP837824PKAZ |
| Abbreviation | matK |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8WJP2 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Prunus laurocerasus</i> (Cherry laurel) (<i>Laurocerasus officinalis</i>) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MEEFQGYLEL DRYQQHDFLH PLIFREYIYA LAHDHGLNRS ILLGYDNKSS LLIIKRLISR MYKQNRFIIS ANDSNQKKNL GYNKNLYSQI ISEGFAVIVE IPFSLRLVSS ATEIVKSYNL RSIHSIFPFL EDKFPQLNYV SDVLIPYPIH LEILVQTLRY WVKDASSLHL LRLFLHEYYN WNSLITSNNF FFFSKSNPRL FLLLYNSHVC EYEFILLFLR NQSSHLQLIS SGIFFERIH F YEKIKYPVEE VFANDFPAAI LWFFKDPFMH YVRYQGKSIL ASKDTPLLMN KWKYLVNLW QCHSYVWSQP GRIYINKLSK HSLDFLGYS SIRPNLSVVR SQMLENSFIT DNAMKKLDTL VPIIPLIGSL AKVKFCNALG HPISKSTWAD SSDFDIIDRF LRIFRNLSHY YSGSSRKSL YRIKYILRLS CVKTLARKHK STVRTFLKRL GSKLLEEFFT EEEQILSLIF PRASYTLKKF YSGRIWYLDI FCINDLINHE |
| Source | Baculovirus |
| Target Names | matK |
| Protein Names | Recommended name: Maturase K Alternative name(s): Intron maturase |
| Expression Region | 1-500 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.