



Recombinant Human Scavenger receptor cysteine-rich domain-containing group B protein (SRCRB4D)

Product Code	CSB-BP837847HU
Abbreviation	SRCRB4D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WTU2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	LQP TPLPFQELRL VGGPSRCRGR LEVMHGGSWG SVCDDDWDVV DANVVCRLQG CGLALPVRP LAFGQGRGPI LLDNVECRGQ EAALSECGR GWGVHNCFH YEDVAVLCDEF LPTQPPTKMLTSRAPPTTL PNGKSEGSVR LVGGANLCQG RVEILHSGLW GTVCDDDWGL PDAAVVCRL GCGAAMAATT NAFFGYGTGH ILLDNVHCEG GEPRLAACQS LGWGVHNCGH HEDAGALCAG LGPPTLTALP SSATREDWAW QTDPSATGVG PQPSRETALL TTAAWAAGKK SGRLRLVGGP GPCRGRVEVL HAGGWGTVCD DDWDFADARV ACREAGCGPA LGATGLGHFG YGRGPVLLDN VGCAGTEARL SDCFHLGWGQ HNCGHEDAG ALCAGPEELG LQVQQDGSET TRVPTPRPRD GHLRLVNGAH RCEGRVELYL GQRWGTVCDD AWDLRAAGVL CRQLGCGQAL AAPGEAHFGP GRGPILLDNV KCRGEESALL LCISHRWDH NCDHSEDASV LCQPS
Source	Baculovirus
Target Names	SSC4D
Protein Names	Recommended name: Scavenger receptor cysteine-rich domain-containing group B protein Alternative name(s): S4-SRCRB
Expression Region	48-575
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.