



# Recombinant Human Scavenger receptor cysteine-rich domain-containing group B protein (SRCRB4D)

<b>Product Code</b>	CSB-YP837847HU
<b>Abbreviation</b>	SRCRB4D
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WTU2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	LQP TPLPFQELRL VGGPSRCRGR LEVMHGGSWG SVCDDDWDVV DANVVCRLQG CGLALPVRP LAFGQGRGPI LLDNVECRGQ EAALSECGR GWGVHNCFH Y EDVAVLCDEF LPTQPTRKM LTRAPPTTL PNGKSEGSVR LVGGANLCQG RVEILHSLW GTVCDDDWGL PDAAVVCRQL GCGAAMAATT NAFFGYGTGH ILLDNVHCEG GEPRLAACQS LGWGVHNCGH HEDAGALCAG LGPPTLTALP SSATREDWAW QTDPSATGVG PQPSRETTALL TTAAWAAGKK SGRLRLVGGP GPCRGRVEVL HAGGWGTVCDD DWDFADARV ACREAGCGPA LGATGLGHFG YGRGPVLLDN VGCAGTEARL SDCFHLGWGQ HNCGHEDAG ALCAGPEELG LQVQQDGSET TRVPTPRPRD GHLRLVNGAH RCEGRVELYL GQRWGTVCDD AWDLRAAGVL CRQLGCGQAL AAPGEAHFGP GRGPILLDNV KCRGEESALL LCISHRWDAH NCDHSEDASV LCQPS
<b>Source</b>	Yeast
<b>Target Names</b>	SSC4D
<b>Protein Names</b>	Recommended name: Scavenger receptor cysteine-rich domain-containing group B protein Alternative name(s): S4-SRCRB
<b>Expression Region</b>	48-575
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.