



Recombinant Human Splicing factor U2AF 26 kDa subunit (U2AF1L4)

Product Code	CSB-MP837850HU
Abbreviation	U2AF1L4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WU68
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEYLASIFGTEKDKVNCSEFYFKIGVCRHGDRCSRLHNKPTFSQEVFTELQEK YGEIEEMNVCDNLGDHLVGNVYVKFRREEDGERAVAELS NRWFNGQAVHGN VPEVASATSCICGPFPRTSRGSSMGGDPGAGHPRGSILATIPERG TIGVPLITG MAASEALAPLPFTPNRDRCSWQDLSSKPPSLSCPILPRLPGSIM
Source	Mammalian cell
Target Names	U2AF1L4
Protein Names	Recommended name: Splicing factor U2AF 26 kDa subunit Alternative name(s): U2 auxiliary factor 26 U2 small nuclear RNA auxiliary factor 1-like protein 4 U2(RNU2) small nuclear RNA auxiliary factor 1-like protein 3 Short name= U2 sma
Expression Region	1-202
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length of Isoform 2
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.