



Recombinant Human UPF0235 protein C15orf40 (C15orf40)

Product Code	CSB-MP837852HU
Abbreviation	C15orf40
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WUR7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLRLRSGLRH LRATPNTRGS ARLLCAEMPK KAGATTKGKS QSKEPERPLP PLGPVAVDPK GCVTIAIHAK PGSKQNAVTD LTAEAVNVAI AAPPSEGEAN AELCRYLSKV LELRKSDVVL DKGKKSREKV VKLLASTTPE EILEKLLKKEA KKT
Source	Mammalian cell
Target Names	C15orf40
Protein Names	Recommended name: UPF0235 protein C15orf40
Expression Region	1-153
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.