



Recombinant Human PIH1 domain-containing protein 2 (PIH1D2)

Product Code	CSB-MP837863HU
Abbreviation	PIH1D2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WWB5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	METSSKGLLT QVTQFWNLLD DLAQSDPEGY EKFIQQQLKE GKQLCAAPEP QLCLQTRILK PKEKILFINL CQWTRIPAPQ STTHPVPLTV GKPEDTTEIS DAYTVIDVAY NPDVLHAAEK DQVKKNQLIQ MAMKCIEEFK QFTLSHSYHI TKFRIKGSIQ RMKQNLMDGIQ TDSIDLREKM RRELTGQIR SSTMSNPDPHF PQLLLPKDQV SGKAVCLIEE ISSTEIQVEM KMPAYELKIV HDHSEKPLKI ELKVELPGIN SVSLCDLSVS EDDLIEVSE KYRLHLNLPK LIDTEMTTAK FIKEKSTLII TMPLV
Source	Mammalian cell
Target Names	PIH1D2
Protein Names	Recommended name: PIH1 domain-containing protein 2
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.