



Recombinant Human Intelectin-2 (ITLN2)

Product Code	CSB-BP837867HU
Abbreviation	ITLN2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WWU7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AAAS SLEMLSREFE TCAFSFSSLP RSCKEIKERC HSAGDGLYFL RTKNGVVYQT FCDMTSGGGG WTLVASVHEN DMRGKCTVGD RWSSQQGNKA DYPEGDGNWA NYNTFGSAEA ATSDDYKNPG YYDIQAKDLG IWHVPNKSPM QHWRNSALLR YRTNTGFLQR LGHNLFGIYQ KYPVKYRSGK CWNDNGPAIP VYDFGDAKK TASYSPYGQ REFVAGFVQF RVFNNERAAN ALCAGIKVTG CNTEHHCIGG GGFFPQGKPR QCGDFSADFWD DGYGTHVKSS CSREITEAAV LLFYR
Source	Baculovirus
Target Names	ITLN2
Protein Names	Recommended name: Intelectin-2 Alternative name(s): Endothelial lectin HL-2
Expression Region	27-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.