



Recombinant Human DNA-binding death effector domain-containing protein 2 (DEDD2)

Product Code	CSB-YP837874HU
Abbreviation	DEDD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WXF8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MALSGSTPAP CWEEDECLDY YGMLSLHRMF EVVGGQLTEC ELELLAFLLD EAPGAAGGLA RARSGLELLL ELERRGQCDE SNLRLLGQLL RVLARHDLLP HLARKRRRPV SPERYSYGTS SSSKRTEGSC RRRRQSSSSA NSQQGQWETG SPPTKRQRRS RGRPSGGARR RRRGAPAAPQ QQSEPARPSS EGKVTCDIRL RVRAEYCEHG PALEQQVASR RPQALARQLD VFGQATAVLR SRDLGSVVD IKFSELSYLD AFWGDYLSGA LLQALRGVFL TEALREAVGR EAVRLLVSVD EADYEAGR RR LLMEEEGGR RPTEAS
Source	Yeast
Target Names	DEDD2
Protein Names	Recommended name: DNA-binding death effector domain-containing protein 2 Alternative name(s): DED-containing protein FLAME-3 FADD-like anti-apoptotic molecule 3
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.