



Recombinant tRNA (Ile)-lysidine synthase (tilS)

Product Code	CSB-BP837932EOD
Abbreviation	tilS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X8W8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	MTLTLNRQLL TSRQILVAFS GGLDSTVLLH QLVQWRTEPN GGALRAIHVH HGLSANADAW VTHCENVCQQ WQVPLVVERV QLAQEGLGIE AQARQARYQA FARTLLPGEV LVTAQHLDDQ CETFLLALKR GSGPAGLSAM AEVSEFAGTR LIRPLLARTR GELAQWALAH GLRWIEDESN QDDSYDRNFL RLRVVPLLQQ RWPFAETTA RSAALCAEQE SLLDELLADD LAHCQSPQGT LQIVPMLAMS DARRAAIIRR WLAGQNAPMP SRDALVRIWQ EVALAREDAS PCLRLGAFEI RRYQSQLWWI KSVTGQSENI VPWQTWLQPL ELPAGQGSVQ LNAGGDIRPP RADEAVSVRF KAPGLLHIVG RGGGRKLKKI WQELGVPPWL RDTTPLLFYG ETLIAAAGGF VTQEGVAEGE NGISFVWQRN A
Source	Baculovirus
Target Names	tilS
Protein Names	Recommended name: tRNA(Ile)-lysidine synthase EC= 6.3.4.19 Alternative name(s): tRNA(Ile)-2-lysyl-cytidine synthase tRNA(Ile)-lysidine synthetase
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.