



Recombinant UDP-N-acetylmuramoyl-L-alanyl-D-glutamate--2,6-diaminopimelate ligase (murE)

Product Code	CSB-EP837940EOD-B
Abbreviation	murE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X9Z2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	ADRNLRDLL APWVPDAPSR ALREMTLDSR VAAAGDLFVA VVGHQADGRR YIPQAIQGV AAIIEAKGE ATDGEIREMH GVPVIYLSQL NERLSALAGR FYHEPSDNLRLVGVGTGTNGK TTTTQLLAQW SQLLGETSAV MGTVGNLLG KVIPTENTTG SAVDVQHELA GLVDQDATFC AMEVSSHGLV QHRVAALKFA ASVFTNLSRD HLDYHGDMEH YEAAKWLLYS EHHCGQAIIN ADDEVGRRWL AKLPDAVAVS MEDHINPNCH GRWLKATEVN YHDSGATIRF SSSWGDGEIE SHLMGAFNVS NLLLALATLL ALGYPLADLL KTAARLQPVC GRMEVFTAPG KPTVVVDYAH TPDALAKALQ AARLHCAGKL WCVFGCGGDR DKGKRPLMGA IAEEFADVAV VTDDNPRTEE PRAIINDILA GMLDAGHAKV MEDRAEAVTC AVMQAKENDV VLVAGKGHED YQIVGNQRLD YSDRVTVARL LGVIA
Source	E.coli
Target Names	murE
Protein Names	Recommended name: UDP-N-acetylmuramoyl-L-alanyl-D-glutamate--2,6-diaminopimelate ligase EC= 6.3.2.13 Alternative name(s): Meso-A2pm-adding enzyme Meso-diaminopimelate-adding enzyme UDP-MurNAc-L-Ala-D-Glu:meso-diaminopimelate ligase
Expression Region	2-495
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.