



Recombinant Chicken Aminomethyltransferase, mitochondrial (AMT)

Product Code	CSB-BP001684CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P28337
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	SAPE GLKQTPLDAL HRARGGRMVP FAGWSLPVQY GRGHLESHLH TRRHCSLFDV SHMLQTRVYG RDRVRFLESL VVGDI AELRP GQGTLLTLTN ERGDIVDDLI VTNTAEDHLY VVSNAGCADK DRAVMEGRAA ELRAAGGDVH LEVSGQRAAG VQGPSMAQVL QAGLPDDLTK LTFMTSTATT VFGVPGCRVT RCGYTGEDGV EISVPAGRAV ELAERLLGCP EVWPAGLAAR DSLRLEAGLC LYGNDIDEST TPVEAGLLWT LGKRRRTAMD FPGAAIIMEQ VKEKPKRKR V GLTSVGPPLR PPAAILGPEG TPVGTVTSGC PSPSLGKNIA MGYVQAAHSR PGTTLTVEVR KKQHPALVTK MPFVPTHYYM AK
Source	Baculovirus
Target Names	AMT
Protein Names	Recommended name: Aminomethyltransferase, mitochondrial EC= 2.1.2.10 Alternative name(s): Glycine cleavage system T protein Short name= GCVT
Expression Region	17-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.