



Recombinant Paralichthys olivaceus Pro-neuropeptide Y (npy)

Product Code	CSB-YP838575ERT
Abbreviation	npy
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q90WF4
Product Type	Recombinant Protein
Immunogen Species	Paralichthys olivaceus (Bastard halibut) (Hippoglossus olivaceus)
Purity	>85% (SDS-PAGE)
Sequence	YP VKPENPGDDA PAEELAKYYS ALRHYINLIT RQRY
Source	Yeast
Target Names	npy
Protein Names	Recommended name: Pro-neuropeptide Y Cleaved into the following 2 chains: 1. Neuropeptide Y Alternative name(s): Neuropeptide tyrosine Short name= NPY C-flanking peptide of NPY Short name= CPON
Expression Region	29-64
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.