



# Recombinant Chicken Protein Wnt-1 (WNT1)

<b>Product Code</b>	CSB-BP838596CH
<b>Abbreviation</b>	WNT1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q91029
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	WWGIVNVASS TLLTDSGSV QLVLEPGLQL LSRKQRRLLIR QNPGILHSVS AGLQSAVREC QWQFRNRRWN CPTSQGPNIIF GKIVNRGCRE TAFIFAITS GVTHSVARSC SEGSIESCTC DYRRRGGPGG DWHWGGCSDN IDFGRLFGRE FVDSSEKGRD LRFLMNLHNN EAGRMTVFSE MRQECKCHGM SGSRTVRTCW MRLPTFRAVG DVLKDRFDGA SRVIYGNKGS NRASRVELHH LEPENPAHKP PSPHDLVYFE KSPNFCTYSG KMGTAGTAGR ACNSSSPGLD GCELLCCGRG FRTRTQRVTE RCNCTFWCC HVSLNCTNT
<b>Source</b>	Baculovirus
<b>Target Names</b>	WNT1
<b>Protein Names</b>	Recommended name: Protein Wnt-1
<b>Expression Region</b>	1-330
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It is very conserved in evolution, and This protein is known to be 98% identical to the mouse Wnt1 protein at the amino acid level. The studies in mouse indicate that the Wnt1 protein functions in the induction of the mesencephalon and cerebellum. This gene was originally considered as a candidate gene for Joubert syndrome, an autosomal recessive disorder with cerebellar hypoplasia as a leading feature. However, further studies suggested that the gene mutations might not have a significant role in Joubert syndrome. This gene is clustered with another family member, WNT10B, in the chromosome 12q13 region.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.