



Recombinant Gloydius ussuriensis Thrombin-like enzyme calobin-1

Product Code	CSB-YP838597GGS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91053
Product Type	Recombinant Protein
Immunogen Species	Gloydius ussuriensis (Ussuri mamushi) (Gloydius blomhoffii ussuriensis)
Purity	>85% (SDS-PAGE)
Sequence	VIGGDE CNINEHRFLV ALYNSRSRTL FCGGTLINQE WVLTAAH CER KNFRIKLG I H SKKVPNEDEQ TRVPKEKFFC LSSKNYTLWD KDIMLIRLDS PVSNSEHIAP LSLPSSPPSV GSVCRIMGWG RISPTKETYP DVPHCANINL LEYEMCRAPY PEFGLPATSR TLCAGILEGG KDTCRGDSGG PLICNGQFQG IASWGDDPCA QPHKPAAYTK VFDHLDWIQS I IAGNTDASC PP
Source	Yeast
Protein Names	Recommended name: Thrombin-like enzyme calobin-1 Short name= SVTLE EC= 3.4.21.- Alternative name(s): Calobin I Fibrinogen-clotting enzyme Snake venom serine protease Short name= SVSP
Expression Region	25-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.