



Recombinant *Xenopus laevis* Fibroblast growth factor 9 (fgf9)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP838635XBE |
| Abbreviation | fgf9 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q91875 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Xenopus laevis</i> (African clawed frog) |
| Purity | >85% (SDS-PAGE) |
| Sequence | LGEVGNV FGVQDAVSFG NVPVLQVDTP VLLSDHMSHH SEAGGLPRGS AVTDLEHLKG ILRRRQLYCR TGFHLEIFPN GTIQGTRQDH NFRGILEFIS IAVGLVSIRG VDSGLYLGMN EKGELYGSEK LTQECVFREQ FEENWYNTYS SNIYKHADTG RRYVALNKD GTSRDGTRTK RHQKFTHFLP RPVDPEKVPE LYKDILSQS |
| Source | E.coli |
| Target Names | fgf9 |
| Protein Names | Recommended name: Fibroblast growth factor 9 Short name= FGF-9 Alternative name(s): Glia-activating factor Short name= GAF HBGF-9 Heparin-binding growth factor 9 XFGF-9 |
| Expression Region | 4-209 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |