



Recombinant Maize streak virus genotype E Replication-associated protein A (C1)

Product Code	CSB-BP838665MIAM
Abbreviation	C1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91MG1
Product Type	Recombinant Protein
Immunogen Species	Maize streak virus genotype E (isolate Pat) (MSV)
Purity	>85% (SDS-PAGE)
Sequence	MASSSSNRSF LHRNANTFLT YPHCPENPEI ISQKLWDLVA RWNPLYIVCA REHRDGNMH LHALLQTDKP VRTTDARIFD IEGFHPNIQS AKSVNKVRDY ILKEPLAVFE RGTIFPRKSC FQGNTPPFPK KNPKNDEIMA HIISHATSKQ EYLCLVRKEF PYDWATKLQY FEYSANKLFP DIQEEFISPH PPSSPDLLCN ESIKDWLQPN IYQVSPEAYM LLQPTCYTVD EAISDLTWMD NLSSQQMKDQ ESRASTSSVQ QGQGNLLGPE V
Source	Baculovirus
Target Names	C1
Protein Names	Recommended name: Replication-associated protein A Short name= RepA EC= 3.1.21.-
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.