



Recombinant Human Retrograde Golgi transport protein RGP1 homolog (RGP1)

Product Code	CSB-MP838781HU
Abbreviation	RGP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92546
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MIEVVAELSR GPVFLAGEAL ECVVTVTNPL PPTATSASSE ALAWASAQIH CQFHASESRV ALPPPDSSQP DVQPDSQTVF LPHRGERGQC ILSTPPKILF CDLRLDPGES KSYSYSEVLP IEGPPSFRGQ SVKYVYKLT I GCQRVNSPIT LLRVPLRVLV LTGLQDVRFP QDEAVAPSSP FLEEDEGGKK DSWLAELAGE RLMAATSCRS LHLYNISDGR GKVGTFGIFK SVYRLGEDVV GTLNLGEGTV ACLQFSVSLQ TEERVQPEYQ RRRGAGGVPS VSHVTHARHQ ESCLHTTTRTS FSLPIPLSST PGFCTAIVSL KWRLHFEFVT SREPLVLLP PVEQPEPTTW TGPEQVPVDT FSWDLPIKVL PTSPTLASYA APGPSTSTIT I
Source	Mammalian cell
Target Names	RGP1
Protein Names	Recommended name: Retrograde Golgi transport protein RGP1 homolog
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.