



Recombinant Human AP-3 complex subunit sigma-1 (AP3S1)

Product Code	CSB-MP838783HU
Abbreviation	AP3S1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92572
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MIKAILIFNN HGKPRLSKFY QPYSED TQQQ IIRETFHLVS KRDENVCNFL EGLLIGGSD NKLIYRHYAT LYFVFCVDSS ESELGILDLI QV FVETLDKC FENVCELDLI FHVDKVHNIL AEMVMGGMVL ETNMNEIVTQ IDAQNKLEKS EAGLAGAPAR AVSAVKNMNL PEIPRNINIG DISIKVPNLP SFK
Source	Mammalian cell
Target Names	AP3S1
Protein Names	Recommended name: AP-3 complex subunit sigma-1 Alternative name(s): AP-3 complex subunit sigma-3A Adapter-related protein complex 3 sigma-1 subunit Clathrin-associated/assembly/adaptor protein, small 3 Sigma-3A-adaptin Short name=
Expression Region	1-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.