



Recombinant Human General transcription factor IIH subunit 4 (GTF2H4)

Product Code	CSB-BP838803HU
Abbreviation	GTF2H4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92759
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MESTPSRGLN RVHLQCRNLQ EFLGGLSPGV LDRLYGHPAT CLAVFRELP LAKNWVMRML FLEQPLPQAA VALWVKKEFS KAQEESTGLL SGLRIWHTQL LPGGLQGLIL NPIFRQNLRI ALLGGGKAWS DDTSQLGPDH HARDVPSLDK YAEERWEVVL HFMVGSPSAA VSQDLAQLLS QAGLMKSTEP GEPPCITSAG FQFLLLDTPA QLWYFMLQYL QTAQSRGMDL VEILSFLFQL SFSTLGKDYS VEGMSDLLN FLQHLREFGL VFQRKRKSRR YYPTRLAINL SSGVSGAGGT VHQPFGIVVE TNYRLYAYTE SELQIALIAL FSEMLYRFPN MVVAQVTRES VQQAIASGIT AQQIIHFLRT RAHPVMLKQT PVLPPITDQ IRLWELERDR LRFTEGVLYN QFLSQVDFEL LLAHARELGV LVFENSAKRL MVVTPAGHSD VKRFRWKRQKH SS
Source	Baculovirus
Target Names	GTF2H4
Protein Names	Recommended name: General transcription factor IIH subunit 4 Alternative name(s): Basic transcription factor 2 52 kDa subunit Short name= BTF2 p52 General transcription factor IIH polypeptide 4 TFIIH basal transcription factor compl
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.