



Recombinant Human Caspase-10 (CASP10)

Product Code	CSB-EP838815HU-B
Abbreviation	CASP10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92851
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	V KTFLEALPQE SWQNKHAGSN GNRATNGAPS LVSRGMQGAS ANTLNSETST KRAAVYRMNR NHRGLCVIVN NHSFTSLKDR QGTHKDAEIL SHVFQWLGFT VHIHNNVTKV EMEMVLQKQK CNPAHADGDC FVFCILTHGR FGAVYSSDEA LIPIREIMSH FTALQCPRLA EKPKLFFIQA CQGEEIQPSV SIEAD
Source	E.coli
Target Names	CASP10
Protein Names	Recommended name: Caspase-10 Short name= CASP-10 EC= 3.4.22.63 Alternative name(s): Apoptotic protease Mch-4 FAS-associated death domain protein interleukin-1B-converting enzyme 2 Short name= FLICE2 ICE- like apoptotic
Expression Region	220-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a protein which is a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein cleaves and activates caspases 3 and 7, and the protein itself is processed by caspase 8. Mutations in this gene are associated with apoptosis defects seen in type II autoimmune lymphoproliferative syndrome. Three alternatively spliced transcript variants encoding different isoforms have been described for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.