



# Recombinant Human Putative deoxyribonuclease TATDN2 (TATDN2)

<b>Product Code</b>	CSB-BP839001HU
<b>Abbreviation</b>	TATDN2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q93075
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MASERGKVKH NWSSTSEGCP RKRSCLEPC DVAPSSRPAQ RSASRSGGPS SPKRLKAQKE DDVACSRRLS WGSSRRRNNS SSSFSPHFLG PGVGGGAASKG CLIRNTRGFL SSGGSPLRPA NASLEEMASL EEEACSLKVD SKDSSHNSTN SEFAAEAEQG NDTIEEPNKV QKRKRDRLRD QGSTMILKA IQGILGKSMP KRKGAAATRA KPSAAEHPSH GEGPARSEGP AKTAEGAARS VTVTAAQKEK DATPEVSMEE DKTVPERSSF YDRRVVIDPQ EKPSEEPDGD RRTVIDKCSP PLEFLDDSDS HLEIQKHKDR EVVMEHPSSG SDWSDVEEIS TVRFSQEEPV SLKPSAVPEP SSFTTDYVMY PPHLYSSPWC DYASYWTSSP KPSSYPSTGS SSNDAAQVGK SSRSRMSDYS PNSTGSVQNT SRDMEASEEG WSQNSRSFRF SRSSEEREVK EKRTFQEEMP PRPCGGHASS SLPKSHLEPS LEEGFIDTHC HLDMLYSKLS FQGTFTKFRK IYSSSFPKEF QGCISDFCDP RLTLDCLWEE LLKEDLVWGA FGCHPHFARY YSESQERNLL QALRHPKAVA FGEMGLDYSY KCTTPVPEQH KVFERQLQLA VSLKKPLVIH CREADEDLLE IMKKFVPPDY KIHRCFTGS YPVIEPLLKY FPNMSVGFTA VLTYSSEAWEA REALRQIPLE RIIVETDAPY FLPRQVPKSL CQYAHPLGLAL HTVREIARVK DQPLSLTLAA LRENTSRLYS L
<b>Source</b>	Baculovirus
<b>Target Names</b>	TATDN2
<b>Protein Names</b>	Recommended name: Putative deoxyribonuclease TATDN2 EC= 3.1.21.-
<b>Expression Region</b>	1-761
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.