



Recombinant Arabidopsis thaliana Germination-specific cysteine protease 1 (GCP1)

Product Code	CSB-MP839115DOA
Abbreviation	GCP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94B08
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	VPETVD WRQKGAVNPI KDQGTGCGSCW AFSTTAAVEG INKIVTGELI SLSEQELVDC DKSYNQGCNG GLMDYAFQFI MKNGLNTEK DYPYRGFGGK CNSFLKNSRV VSIDGYEDVP TKDETALKKA ISYQPVSAI EAGGRIFQHY QSGIFTGSCG TNLDAHVVAV GYGSENGVDY WIVRNSWGPR WGEEGYIRME RNLAASKSGK CGIAVEASYP VKYSPNPVVRG NTISSV
Source	Mammalian cell
Target Names	GCP1
Protein Names	Recommended name: Germination-specific cysteine protease 1 EC= 3.4.22.-
Expression Region	145-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.