



Recombinant Transcription factor che-1 (che-1)

Product Code	CSB-BP839268CXY
Abbreviation	che-1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q966L8
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSYHYQMNNN NTDFFDSDANY IPQMQMHHGY QDFYASDSTS SMPSSGTTQQ ALQQTENGII QGQNGSDHGV AFHFRAYGGW RQHDWNSTSS WSIGNQLEME NGTLTMSHGA TNDGNIRNNQ QMMTSNHITT DSFRQIPPSG LCYMDPTPSL RLPKKEPLPV QRPMARSTPK PFRCQTCGKA FSQAANLTAH KRIHTGEKPF MCPVCNRPFS QSSSLVTHRR THTGERPYPC AQCEKAFTDS STLTKHLRTH TGHKPYVCSI CMMKFTQSGN LHRHMKTHK
Source	Baculovirus
Target Names	che-1
Protein Names	Recommended name: Transcription factor che-1 Alternative name(s): Abnormal chemotaxis protein 1
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>