



Recombinant Human Lethal (3)malignant brain tumor-like protein 2 (L3MBTL2)

Product Code	CSB-MP839279HU
Abbreviation	L3MBTL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q969R5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEKPRSIEET PSSEPMEEEE DDDLELFGGY DSFRSYNSSV GSESSSYLEE SSEAENEDRE AGELPTSPLH LLSPGTPrSL DGSGSEPAVC EMCGIVGTRE AFFSKTKRFC SVSCRSYSS NSKKASILAR LQGKPPTKKA KVLHKAAWSA KIGAF LHSQG TGQLADGTPT GQDALVLGFD WGKFLK DHSY KAAPVSCFKH VPLYDQWEDV MKGMKVEVLN SDAVLPSRVY WIASVIQTAG YRVLLRYEGF ENDASHDFWC NLGTVDVHPI GWCAINSKIL VPPRTIHAKF TDWKGYLMKR LVGSRTL PVD FHIKMVESMK YPFRQGMRL E VVDKSQVSR T RMAVVDTVIG GRLRLLYEDG DSDDDFWCHM WSPLIHPVGW SRRVGHGIKM SERRSDMAHH PTFRKIYCDA VPYLFKKVRA VYTEGGWFEE GMKLEAIDPL NLGNICVATV CKVLLDGYLM ICVDGGPSTD GLDWFCYHAS SHAI FPATFC QKNDIELTPP KGYEAQTFNW ENYLEKTKSK AAPSR LFNMD CPNHGFKVGM KLEAVDLMEP RLICVATVKR VVHRLLSIHF DGWDSEYDQW VDCESPDIYP VGWCELTGYQ LQPPVAAEPA TPLKAKEATK KKKKQFGKKR KRIPPTKTRP LRQGSKKPLL EDDPQGARKI SSEPVPGEII AVR VKEEHL D VASPKASSP ELPVSVENIK QETDD
Source	Mammalian cell
Target Names	L3MBTL2
Protein Names	Recommended name: Lethal(3)malignant brain tumor-like protein 2 Short name= H-l(3)mbt-like protein 2 Short name= L(3)mbt-like protein 2
Expression Region	1-705
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.