



Recombinant Human Lethal (3)malignant brain tumor-like protein 2 (L3MBTL2)

| | |
|--------------------------|--|
| Product Code | CSB-YP839279HU |
| Abbreviation | L3MBTL2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q969R5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MEKPRSIET PSSEPMEEEE DDDLELFGGY DSFRSYNSSV GSESSSYLEE SSEAENEDRE AGELPTSPLH LLSPGTPrSL DGSGSEPAVC EMCGIVGTRE AFFSKTKRFC SVSCRSYSS NSKKASILAR LQGKPPTKKA KVLHKAWSA KIGAFHLSQG TGQLADGTPT GQDALVLGFD WGKFLKDHSY KAAPVSCFKH VPLYDQWEDV MKGMKVEVLN SDAVLPSRVY WIASVIQTAG YRVLLRYEGF ENDASHDFWC NLGTVDVHPI GWCAINSKIL VPPRTIHAKF TDWKGylMKR LVGSRTLpVD FHIKMVESMK YPFRQGMRLV VDKSQVSRt RMAVVDtVIG GRLRLLYEDG DSDDDFWCHM WSPLIHPVGW SRRVGHGIKM SERRSDMAHH PTFRKIYCDA VPYLFKKVRA VYTEGGWFEE GMKLEAIDPL NLGNICVATV CKVLLDGyLM ICVDGGPSTD GLDWFCYHAS SHAIpATFC QKNdIELTPP KGYEAQTFNW ENYLEKTKSK AAPSRlFNMD CPNHGfKVGm KLEAVDLMEP RLICVATVKR VVHRLLSIHf DGWDSEYDQW VDCESPDIYP VGWCELTGYQ LQPPVAAEPA TPLKAKEATK KKKKQFGKKR KRIPPTKTRP LRQGSKKPLL EDDPQGARKI SSEPVPGEII AVRvKEEHLd VASpDKASSP ELPVSVENIK QETDD |
| Source | Yeast |
| Target Names | L3MBTL2 |
| Protein Names | Recommended name: Lethal(3)malignant brain tumor-like protein 2 Short name= H-l(3)mbt-like protein 2 Short name= L(3)mbt-like protein 2 |
| Expression Region | 1-705 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.