



Recombinant Human RILP-like protein 2 (RILPL2)

Product Code	CSB-MP839282HU
Abbreviation	RILPL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q969X0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEEPPVREEE EEEGEEDEER DEVGPEGALG KSPFQLTAED VYDISYLLGR ELMALGSDPR VTQLQFKVVR VLEMLEALVN EGSLALEELK MERDHLRKEV EGLRRQSPPA SGEVNLGPNK MVVDLTDPNR PRFTLQELRD VLQERNKLKS QLLVVQEELQ CYKSGLIIPR EGPGRREKD AVVTSAKNAG RNKEEKTIK KLFFFRSGKQ T
Source	Mammalian cell
Target Names	RILPL2
Protein Names	Recommended name: RILP-like protein 2 Alternative name(s): Rab-interacting lysosomal protein-like 2 p40phox-binding protein
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.