



Recombinant Human Coiled-coil domain-containing protein 102A (CCDC102A)

Product Code	CSB-EP839284HU-B
Abbreviation	CCDC102A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96A19
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSHGSPRLA ESPQLSKGSL LTILGSPSPE RMGPADSLPP TPPSGTPSPG PPPALPLPPA PALLADGDWE SREELRLREL EEARARAAQM EKTMRWSDC TANWREKWSK VRAERNRARE EVRQLRQRLD ALTKELAGAR RERQEAQGEC EARGRELARL RGARGVADQT RDGPEPEAER EPVRDVGSER PPGSQELELV ESKLKSMPPEE SEDCWEARSL GAGGPRGSSG RQERSRLPWE DTAATEEEAS KLTALRLRLD ESQKVLLKER EDKLALSRNI EKLEGELSQW KIKYEELSKT KQEMLKQLSI LKEAHQDELG RMSEDELEDEL GARSSMDRKM AELRGEMERL QAENAAEWGR RERLETEKLG LERENKKLRA QVGDLEEALA RRRRQTASAL DCDLRASQAA LFEKNKELAD LKHVHGKLLK QFQEKVAELA HANRRVEQHE AEVKKLRLRV EELKKELAQA EDELDEAHNQ ARKLQRSLE QTEQSENLQV QLEHLQSRLR RQQQNAPLFG KIRSARFGTE EAEDGTSDDL EDELQIQVA
Source	E.coli
Target Names	CCDC102A
Protein Names	Recommended name: Coiled-coil domain-containing protein 102A
Expression Region	1-550
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.