



# Recombinant Human Coiled-coil domain-containing protein 102A (CCDC102A)

<b>Product Code</b>	CSB-YP839284HU
<b>Abbreviation</b>	CCDC102A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96A19
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSHGSPRLA ESPQLSKGSL LTILGSPSPE RMGPADSLPP TPPSGTPSPG PPPALPLPPA PALLADGDWE SREELRLREL EEARARAAQM EKTMRRWSDC TANWREKWSK VRAERNRARE EVRQLRQRLD ALTKELAGAR RERQEAQGEC EARGRELARL RGARGVADQT RDGPEPEAER EPVRDVGSER PPGSQELELV ESKLKSMPPEE SEDCWEARSL GAGGPRGSSG RQERSRLPWE DTAATEEEAS KLTALRLRLD ESQKVLLKER EDKLALSRNI EKLEGELSQW KIKYEELSKT KQEMLKQLSI LKEAHQDELG RMSEDELEDEL GARSSMDRKM AELRGEMERL QAENAAEWGR RERLETEKLG LERENKKLRA QVGDLEEALA RRRRQTASAL DCDLRASQAA LFEKNKELAD LKHVHGKLLK QFQEKVAELA HANRRVEQHE AEVKKLRLRV EELKKELAQA EDELDEAHNQ ARKLQRSLE QTEQSENQV QLEHLQSRLR RQQNAPLFG KIRSARFGTE EAEDGTSDLD EDELQIQVA
<b>Source</b>	Yeast
<b>Target Names</b>	CCDC102A
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 102A
<b>Expression Region</b>	1-550
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.