



Recombinant Human Adrenocortical dysplasia protein homolog (ACD)

Product Code	CSB-MP839289HU
Abbreviation	ACD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96AP0
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAGSGRLVLRPWIRELILGSETPSSPRAGQLLEVLQDAEAAVAGPSHAPDTS VGATLLVSDGTHSVRCLVTREALDTSWEEKEFGFRGTEGRLLLLQDCGVHV QVAEGGAPAEFYLVQDRFSLLPTEQPRLRVPGCNQDLVDVQKKLYDCLEEHLS ESTSSNAGLSLSQLLDEMREDQEHQALVCLAESCLTLEGPCTAPPVTHWAA SRCKATGEAVYTPSSMLCISENDQLILSSLGPCQRTQGPELPPDPALQDLS LTLIASPPSSPSSSGTPALPGHMSSEESGTSISLLPALSLAAPDPGQRSSSQPS PAICSAPATLTPRSPHASRTPSSPLQSCTPSLSRSHVPSPHQALVTRPQKPS LEFKEFVGLPCKNRPPFPRTGATRGAQEPCSVWEPPKRHRDGSFAFQYIEYEP PCTSLCARVQAVRLPPQLMAWALHFLMDAQPGSEPTPM
Source	Mammalian cell
Target Names	ACD
Protein Names	Recommended name: Adrenocortical dysplasia protein homolog Alternative name(s): POT1 and TIN2-interacting protein
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.