



Recombinant Human Protein FAM131C (FAM131C)

Product Code	CSB-YP839290HU
Abbreviation	FAM131C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96AQ9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGSCVSRDLF TSAHKNCMP QGADPLNPD LPSGRTP TVAP DCVIGKDKQM DFCWDPWQRC FQTTNGYLS DSRSRPGNY NV AALATSSLVG VVQSIKDHIT KPTAMARGRV AHLIEWKGWS AQPAGWELSP AEDEHYCCLP DELREARFAA GVAEQFAITE ATLSAWSSLD EEELHPENSP QGIVQLQDLE SIYLQDSLPS GPSQDDSLQA FSSPSPSPDS CPSPEEPPST AGIPQPPSPE LQHRRRLPGA QGPEGGTHPP GSLPSMDSGS LWEEDEVFYN
Source	Yeast
Target Names	FAM131C
Protein Names	Recommended name: Protein FAM131C
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.