



Recombinant Human Arrestin domain-containing protein 3 (ARRDC3)

Product Code	CSB-BP839293HU
Abbreviation	ARRDC3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96B67
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVLGKVKSLT ISFDCLNDSN VPVYSSGDTV SGRVNLEVTG EIRVKSLKIH ARGHAKVRWT ESRNAGSNTA YTQNYTEEVE YFNHKDILIG HERDDDNSEE GFHTIHSGRH EYAFSFLPQ TPLATSFEGR HGSVRYWVKA ELHRPWLLPV KLLKKEFTVFE HIDINTPSLL SPQAGTKEKT LCCWFCTSGP ISLSAKIERK GYTPGESIQI FAEIENCSSR MIVPKAAIQ TQAFYAKGKM KEVKQLVANL RGESLSSGKT ETWNGKLLKI PPVSPSILDC SIIRVEYS LM VYVDIPGAMD LFLNLPLVIG TIPLHPFGSR TSSVSSQCSM NMNWLSLSLP ERPEAPPSYA EVTTEEQR RN NLAPVSACDD FERALQGPLF AYIQEFRFLP PPLYSEIDPN PDQSADDRPS CPSR
Source	Baculovirus
Target Names	ARRDC3
Protein Names	Recommended name: Arrestin domain-containing protein 3 Alternative name(s): TBP-2-like inducible membrane protein Short name= TLIMP
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.