



# Recombinant Human Arrestin domain-containing protein 3 (ARRDC3)

<b>Product Code</b>	CSB-YP839293HU
<b>Abbreviation</b>	ARRDC3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96B67
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVLGKVKSLT ISFDCLNDSN VPVYSSGDTV SGRVNLEVTG EIRVKSLKIH ARGHAKVRWT ESRNAGSNTA YTQNYTEEVE YFNHKDILIG HERDDDNSEE GFHTIHSGRH EYAFSFLPQ TPLATSFEGR HGSVRYWVKA ELHRPWLLPV KLLKKEFTVFE HIDINTPSLL SPQAGTKEKT LCCWFCTSGP ISLSAKIERK GYTPGESIQI FAEIENCSSR MVVPKAAIQ TQAFYAKGKM KEVKQLVANL RGESLSSGKT ETWNGKLLKI PPVSPSILDC SIIRVEYS LM VYVDIPGAMD LFLNLPLVIG TIPLHPFGSR TSSVSSQCSM NMNWLSLSLP ERPEAPPSYA EVTTEEQR RN NLAPVSACDD FERALQGPLF AYIQEFRFLP PPLYSEIDPN PDQSADDRPS CPSR
<b>Source</b>	Yeast
<b>Target Names</b>	ARRDC3
<b>Protein Names</b>	Recommended name: Arrestin domain-containing protein 3 Alternative name(s): TBP-2-like inducible membrane protein Short name= TLIMP
<b>Expression Region</b>	1-414
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.