



Recombinant Human Protein FAM149B1 (FAM149B1)

Product Code	CSB-EP839295HU-B
Abbreviation	FAM149B1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BN6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MISRYTRKAV PQSLELKGIT KHALNHHPPP EKLEEISPTS DSHEKDTSSQ SKSDITRESS FTSADTGNSL SAFPSYTGAG ISTEKSSDFS WGYGELDQNA TEKVQTMFTA IDELLYEQKL SVHTKSLQEE CQQWTASFPH LRILGRQIIT PSEGYRLYPR SPSAVSASYE TTLSQERDST IFGIRGKKLH FSSSYAHKAS SIAKSSSFCS MERDEEDSII VSEGIIEEYL AFDHIDIEEG FHGKKSEAAT EKQKLGYPPI APFYCMKEDV LAYVFDVWC KVVSCMEQLT RSHWEGFASD DESNVAVTRP DSESSCVLSE LHPLVLPVP QSKVLYITSN PMSLCQASRH QPNVNDLLVH GMPLQPRNLS LMDKLLDLD KLLMRPGSST ILSTRNWPNR AVEFSTSSLS YTVQSTRRRN PPPRTLHPIS TSHSCAETPR SVEEILRGAR VPVAPDSLSS PSPTPLSRNN LLPIGTAEV EHVSTVGPQR QMKPHGDSSR AQSAVVDEPN YQQPQERLLL PDFFPRPNTT QSFLDQTQYR RSCAVEYPHQ ARPGRGSAGP QLHGSTKSQS GGRPVSRTTRQ GP
Source	E.coli
Target Names	FAM149B1
Protein Names	Recommended name: Protein FAM149B1
Expression Region	1-582
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.