



Recombinant Human Protein FAM149B1 (FAM149B1)

Product Code	CSB-YP839295HU
Abbreviation	FAM149B1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BN6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MISRYTRKAV PQSLELKGIT KHALNHHPPP EKLEEISPTS DSHEKDTSSQ SKSDITRESS FTSADTGNSL SAFPSYTGAG ISTEGLSSDFS WGYGELDQNA TEKVQTMFTA IDELLYEQKL SVHTKSLQEE CQQWTASFPH LRILGRQIIT PSEGYRLYPR SPSAVSASYE TTLSQERDST IFGIRGKKLH FSSSYAHKAS SIAKSSSFCFS MERDEEDSII VSEGIIEEYL AFDHIDIEEG FHGKKSEAAT EKQKLGYPPI APFYCMKEDV LAYVFDVSWC KVVSCMEQLT RSHWEGFASD DESNVAVTRP DSESSCVLSE LHPLVLPVVP QSKVLYITSN PMSLCQASRH QPNVNDLLVH GMPLQPRNLS LMDKLLDLDLDD KLLMRPGSST ILSTRNWPNR AVEFSTSSLS YTVQSTRRRR PPPRTLHPIS TSHSCAETPR SVEEILRGAR VPVAPDSLSS PSPTPLSRNN LLPPIGTAEV EHVSTVGPQR QMKPHGDISSR AQSAVVDEPN YQQPQERLLL PDFFPRPNTT QSFLLDLDTQYR RSCAVEYPHQ ARPGRGSAGP QLHGSTKSQS GGRPVSRTTRQ GP
Source	Yeast
Target Names	FAM149B1
Protein Names	Recommended name: Protein FAM149B1
Expression Region	1-582
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.