



Recombinant Human TRAF-interacting protein with FHA domain-containing protein A (TIFA)

Product Code	CSB-MP839301HU
Abbreviation	TIFA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96CG3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTSFEDADTE ETVTCLQMTV YHPGQLQCGI FQSISFNREK LPSSEVVKFG RNSNICHYTF QDKQVSRVQF SLQLFKKFNS SVLSFEIKNM SKKTNLIVDS RELGYLNKMD LPYRCMVRFG EYQFLMEKED GESLEFFETQ FILSPRSLQ ENNWPPIRPI PEYGTYSLCS SQSSSPTMD ENES
Source	Mammalian cell
Target Names	TIFA
Protein Names	Recommended name: TRAF-interacting protein with FHA domain-containing protein A Alternative name(s): Putative MAPK-activating protein PM14 Putative NF-kappa-B-activating protein 20 TRAF2-binding protein
Expression Region	1-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.