



Recombinant Human Probable RNA polymerase II nuclear localization protein SLC7A6OS (SLC7A6OS)

Product Code	CSB-EP839303HU
Abbreviation	SLC7A6OS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96CW6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEAARTAVLR VKRKRSAEPA EALVLACKRL RSDAVESAAQ KTSEGLERAA ENNVFHLVAT VCSQEEPVQP LLREVL RPSR DSQQRVRRNL RASAREVRQE GRYRVLSSRR SLGTTSSGQE SEYTPGNPEA AGNSGFQLLD LVHEEGEPEA ASAGSCKTSD PDVILCNSVE LIRERLTVSE DPGGVRRQEE QKHDDYVYDI YYLETATPGW IENILSVQPY SQEWELVNDQ QEPEDIYDDE DDENSENNWR NEYPEEESSD GDEDSRGSAD YNSLSEEERG SSRQRMWSKY PLDVQKEFGY DSPHDLDS
Source	E.coli
Target Names	SLC7A6OS
Protein Names	Recommended name: Probable RNA polymerase II nuclear localization protein SLC7A6OS Alternative name(s): ADAMS proteinase-related protein Solute carrier family 7 member 6 opposite strand transcript
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.