



Recombinant Human Protein Hook homolog 2 (HOOK2)

Product Code	CSB-MP839313HU
Abbreviation	HOOK2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96ED9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSVDKAEELCG SLLTWLQTFH VPSPCASPQD LSSGLAVAYV LNQIDPSWFN EAWLQGISD PGPNWKLKVS NLKMLVRLSLV EYSQDVLHP VSEEHLDPVS LIGEFSDPAE LGKLLQLVLG CAISCEKKQD HIQRIMTLEE SVQHVVMEDIA QELMTKDTPD SLSPETYGNF DSQSRRYYFL SEEAEEGDEL QQRCLDLERQ LMLLSEEKQS LAQENAGLRE RMGRPEGEGT PGLTAKKLLL LQSQLEQLQE ENFRLESGRE DERLRCAELE REVAELQHRN QALTSQAQA QALKDEMDEL RQSSERAGQL EATLTSCRRR LGELRELRRQ VRQLEERNAG HAERTRQLED ELRRAGSLRA QLEAQRQVQ ELQGQRQEEA MKAEKWLFEC RNLEEKYESV TKEKERLLAE RDSLREANEE LRCAQLQPRG LTQADPSLDP TSTPVDNLAA EILPAELRET LLRLQLENKR LCRQEAADRE RQEELQRHLE DANRARHGLE TQHRLNQQQL SELRAQVEDL QKALQEQGGK TEDAISILLK RKLEEHLQKL HEADLELQRK REYIEELEPP TDSSTARRIE ELQHNLQKKD ADLRAMEERY RRYVDKARMV MQTMEPKQRP AAGAPPELHS LRTQLRERDV RIRHLEMDFE KRSRQREQEE KLLISAWYNM GMALQQRAGE ERAPAHAQSF LAQQLATNS RRGPLGRLAS LNLRPDTKH
Source	Mammalian cell
Target Names	HOOK2
Protein Names	Recommended name: Protein Hook homolog 2 Short name= h-hook2 Short name= hHK2
Expression Region	1-719
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.