



Recombinant Human Protein KTI12 homolog (KTI12)

Product Code	CSB-MP839314HU
Abbreviation	KTI12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96EK9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPLVVF CGLP YSGKSRRAEE LRVALAAEGR AVYVVDDAAV LGAEDPAVYG DSAREKALRG ALRASVERRL SRHDVVILDS LNYIKGFRIE LYCLARAART PLCLVYCVRP GGPIAGPQVA GANENPGRNV SVSWRPRAEE DGRAQAAGSS VLRELHTADS VVNGSAQADV PKELEREESG AAESPALVTP DSEKSAKHGS GAFYSPELLE ALTLRFEAPD SRNRWDRPLF TLVGLLEPLP LAGIRSALFE NRAPPPHQST QSQPLASGSF LHQLDQVTSQ VLAGLMEAQK SAVPGDLLTL PGTTEHLRFT RPLTMAELSR LRRQFISYTK MHPNNENLPQ LANMFLQYLS QSLH
Source	Mammalian cell
Target Names	KTI12
Protein Names	Recommended name: Protein KTI12 homolog
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.