



# Recombinant Human G antigen family D member 3 (XAGE2)

<b>Product Code</b>	CSB-BP839331HU
<b>Abbreviation</b>	XAGE2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96GT9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSWRGRSTYR PRPRRSLQPP ELIGAMLEPT DEEPKEEKPP TKSRNPTPDQ KREDDQGAAE IQVPDLEADL QELCQTKTGD GCEGGTDVKG KILPKAEHFK MPEAGEGKSQ V
<b>Source</b>	Baculovirus
<b>Target Names</b>	XAGE2
<b>Protein Names</b>	Recommended name: G antigen family D member 3 Alternative name(s): Cancer/testis antigen 12.2 Short name= CT12.2 Protein XAGE-2
<b>Expression Region</b>	1-111
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene is a member of the XAGE subfamily, which belongs to the GAGE family. The GAGE genes are expressed in a variety of tumors and in some fetal and reproductive tissues. This gene is strongly expressed in normal testis, and in Ewing s sarcoma, rhabdomyosarcoma, a breast cancer and a germ cell tumor. This protein shares a sequence similarity with other GAGE/PAGE proteins. Because of the expression pattern and the sequence similarity, this protein also belongs to a family of CT (cancer-testis) antigens.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.