



Recombinant Human Oxidoreductase NAD-binding domain-containing protein 1 (OXNAD1)

Product Code	CSB-YP839337HU
Abbreviation	OXNAD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96HP4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AIR IEAASLRLTL STLRHLTLTS IMKSKRKT DH MERTASVLRR EIVSAAKVCG AAESPSVKSL RLLVADQDF SFKAGQWVDF FIPGVS VVGG FSICSSPRLL EQERVIELAV KYTNHPPALW VHNTCTLDCE VAVRVGGEFF FDPQPADASR NLVLIAGGVG INPLLSILRH AADLLREQAN KRNGYEIGTI KLFYSAKNTS ELLFKKNILD LVNEFPEKIA CSLHVTKQTT QINAELKPYI TEGRITEKEI RDHISKETLF YICGPPMTD FFSKQLENNH VPKEHICFEK WW
Source	Yeast
Target Names	OXNAD1
Protein Names	Recommended name: Oxidoreductase NAD-binding domain-containing protein 1 EC= 1.-.-.-
Expression Region	18-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.