



Recombinant Human Uncharacterized protein C7orf34 (C7orf34)

Product Code	CSB-MP839358HU
Abbreviation	C7orf34
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96L11
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SP GPKDGSQTEK TPSADQNQEQ FEEHFVASSV GEMWQVVDMA QQEEDQSSKT AAVHKHSFHL SFCFSLASVM VFSGGPLRRT FPNIQLCFML TH
Source	Mammalian cell
Target Names	LLCFC1
Protein Names	Recommended name: Uncharacterized protein C7orf34 Alternative name(s): MSSP-binding protein CTM-1
Expression Region	29-122
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.