



Recombinant Human Testis-expressed sequence 34 protein (TEX34)

Product Code	CSB-MP839363HU
Abbreviation	TEX34
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96LK8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGVTGAHGFP CCGKGSVEVA EMRDDLSQHQ IQEEQELEAD MLEQKPQLQV DLDLDPDPDP DPELEIGQVP ALLESELYPA LKLEAELDTE ANSNEESDFE EPMQLVCKIE SVHSNMGLPT PQTFRPWSLN SNCRSFTEEN HVSACHHSIS AQTSKHLFWA NKLIQASEHS LQRAINMQLN NGSAGQPIRS PLREAIPTNA LCSEEQLQIP DAHSAPPTTS SQAPSPLLSS DLPPPIDLTE LITFASSLAM ASSSRMDLPS LEHMMKAPPQ EALEPSTEPL LTTVEEREPE NHAETLPEKP REARAPLKSW SQEDKNFAQS YFDFSKPGIK RATIKGQIQL LQPPATSPLL QGSKEDSVPP GKEKENPLLK KIHFKLSAPT IPEK
Source	Mammalian cell
Target Names	SPATA32
Protein Names	Recommended name: Testis-expressed sequence 34 protein
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.