



Recombinant Human Coiled-coil domain-containing protein 7 (CCDC7)

Product Code	CSB-EP839370HU
Abbreviation	CCDC7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96M83
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKPVKHLTT SNKSANVPAL TTKKGLHNLP LSPELKEKHN AKLIHDKIEP MVLRSPTGE SILRYALPIP SSKTKNLLPE DEMIGKIIKH LKMVVSTLEE TYGHCDQNGE EPFVKHEHEE LLSVGDMDN SFLTYCSQFA AQLLEEALKEE QNILESLFKW FQWQVNQMEE ISKDQTLQA EPPKPKDQVI LNIAEIVRLV QRFEELKNRL KQRSSSVKV MLSKTMDKEN RPEAVKSCEA LAQKIEEFLE AHSTDEFKDV SATEPQTAHS MTNRFNAMLK VFENQANMLE RAVNDQVLLD AEYKQMOCDF QLLSEEKLV L ENELQKLKDK ETKTPTNRT KKAVKTVKKK DKGKSEDSEK KMSPEKEFKI KEDLDQVQKV ARLEIENKVL QEQLKQALQE AEKAKHQLNY FLNQEKLLKS EGKTETTMQV GNSQTKVKGE DSKNIPLEKE TRKSLVSDSG GQRTSDKIQE YPQITAQSGR LIEKRC
Source	E.coli
Target Names	CCDC7
Protein Names	Recommended name: Coiled-coil domain-containing protein 7
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.