



Recombinant Human SH3 and cysteine-rich domain-containing protein 3 (STAC3)

Product Code	CSB-MP839373HU
Abbreviation	STAC3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96MF2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTEKEVLESP KPSFPAETRQ SGLQRLKQLL RKGSTGTKEM ELPPEPQANG EAVGAGGGPI YYIYEEEEEE EEEEEPPPE PPKLVNDKPH KFKDHFFKKP KFCDVCARMI VLNNKFGLRC KNCKTNIHEH CQSYVEMQRC FGKIPPGFHR AYSSPLYSNQ QYACVKDLSA ANRNDPVFET LRTGVIMANK ERKKGQADKK NPVAAMMEEE PESARPEEGK PQDGNPEGDK KAEKKTTPDDK HKQPGFQQSH YFVALYRFKA LEKDDLDLDFPP GEKITVIDDS NEEWWRGKIG EKVGFPPNF IIRVRAGERV HRVTRSFVGN REIGQITLKK DQIVVQKGDG AGGYVKVYTG RKVGLFPTDF LEEI
Source	Mammalian cell
Target Names	STAC3
Protein Names	Recommended name: SH3 and cysteine-rich domain-containing protein 3
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.