



Recombinant Human Protein FAM161B

(FAM161B)

Product Code	CSB-MP839378HU
Abbreviation	FAM161B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96MY7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTVGRPEGAP GGAEGSRQIF PPESFADTEA GEELSGDGLV LPRASKLDEF LSPEEEIDST SDSTGSIYQN LQELKQKGRW CLLES LFQSD PESDENLSED EEDLESFFQD KDRGMVQVQC PQALRCGSTR RCSSLNNLPS NIPRPQTQPP SGSRPPSQHR SVSSWASSIT VPRPFRMTLR EARKKAEWLG SPASFEQERQ RAQRQGEEEA ECHRQFRAQP VPAHVYLYLY QEIMERSEAR RQAGIQKRKE LLLSSLKPFS FLEKEEQLKE AARQRDLAAT AEAKISKQKA TRIPKSI PALGDKLQEA ELFRKIRIQM RALDMLQMAS SPIASSSNRA NPQPRTATRT QQEKLGLFHT NFRFQPRVNP VVPDYEGLYK AFQRRRAAKRR ETQEATRNKP FLLRTANLRH PQRPCDAATT GRRQDSPQPP ATPLPRSRSL SGLASLSANT LPVHITDATR KRESAVRSAL EKKNKADESI QWLEIHKKKS QAMSKSVTLR AKAMDPHKSL EEVFKAKLKE NRNNDKRKRAK EYKKELEEMK QRIQTRPYLF EQVAKDLAKK EAEQWYLDLTL KQAGLEEDFV RNKGQGTRAV QEKETKIKDF PRFQETTKLS IRDPEQGLEG SLEQPASPRK VLEELSHQSP ENLVSLA
Source	Mammalian cell
Target Names	FAM161B
Protein Names	Recommended name: Protein FAM161B
Expression Region	1-647
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.