



Recombinant Human Uncharacterized protein C22orf13 (C22orf13)

Product Code	CSB-MP839384HU
Abbreviation	C22orf13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96NT3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRTEAEAAGP PLEPGDFVQL PVPVIQQLYH WDCGLACSRM VLRYLGLQDD SEFERALQKL QLTRSIWTID LAYLMHHFGV RHRFCTQTLG VDKGYKNQSF YRKHFDTEET RVNQLFAQAK ACKVLVEKCT VSVKDIQAH L A QGHVAIVLV NSGVLHCDLC SSPVKYCCFT PSGHHCF CRT PDYQGHFIVL RGYNRATGCI FYNNPAYADP GMCSTISISNF EEARTSYGTD EDILFVYLD S
Source	Mammalian cell
Target Names	GUCD1
Protein Names	Recommended name: Uncharacterized protein C22orf13 Alternative name(s): Protein LLN4
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.