



Recombinant Human NudC domain-containing protein 1 (NUDCD1)

Product Code	CSB-MP839409HU
Abbreviation	NUDCD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96RS6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEVAANCSLR VKRPLLDPRF EGYKLSLEPL PCYQLELDAA VAEVKLRDDQ YTLEHMHAFG MYNYLHCDSW YQDSVYYIDT LGRIMNLTVM LDTALGKPRE VFRLPTDLTA CDNRLCASIHFSSSTWVTLSDGTGRLYVIG TGERGNSASE KWEIMFNEEL GDPFIIHSI SLLNAEEHSI ATLLLRIEKE ELDMMKSGSFY VSLEWVTISK KNQDNKKYEI IKRDILRGKS VPHYAAIEPD GNGLMIVSYK SLTFVQAGQD LEENMDEDIS EKIKEPLYYW QQTEDDLTVT IRLPEDSTKE DIQIQLPDH INIVLKDHFQ LEGKLYSSID HESSTWIIKE SNSLEISLIK KNEGLTWPEL VIGDKQGELI RDSAQCAAIA ERLMHLTSEE LNPNDKEKP PCNAQELEEC DIFFEESSSL CRFDGNTLKT THVVNLGSNQ YLFSVIVDPK EMPCFCLRHD VDALLWQPHS SKQDDMWEHI ATFNALGYVQ ASKRDKKFFA CAPNYSYAAL CECLRRVFIY RQPAPMSTVL YNRKEGRQVG QVAKQQVASL ETNDPILGFQ ATNERLFVLT TKNLFLIKVN TEN
Source	Mammalian cell
Target Names	NUDCD1
Protein Names	Recommended name: NudC domain-containing protein 1 Alternative name(s): Chronic myelogenous leukemia tumor antigen 66 Tumor antigen CML66
Expression Region	1-583
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.